



**Overview**

Calcium in lactose modified milk powder is determined by KAP analysis (multiple known addition) on 960 (930) Meter. Aliquots of the calcium standard are added automatically to a sample in water. A calcium electrode is used to measure changes in potential after each addition and calculate the sample concentration.

<b>Market</b>	Food and Beverage	<b>Species Measured</b>	Calcium
<b>Sample</b>	Milk powder, lactose modified	<b>Sample Size</b>	0.025 g
		<b>Typical Concentration</b>	5000 mg/100g
<b>Technique #</b>	2 Multiple Known Addition	<b>Electrode</b>	Ca 9320BN; SJ Ref 900100; pH electrode, e.g., 8172BN or 9172BN

**Solutions** 0.1 M Calcium standard 922006; Calcium ISA 932011; Reference electrode filling solution 900011; Concentrated Nitric acid (15.8 M); pH buffer 4.01 910104 and pH buffer 7.00 910107.

**Sample Prep** Step 1: Weigh about 0.25 g sample (record exact weight of sample) in a weigh dish. Transfer to a 500 mL volumetric flask partly filled with 0.001 M HNO<sub>3</sub>. Add 10 mL of the calcium ISA. Fill to the mark with the 0.001 M HNO<sub>3</sub> and mix.  
Step 2: Pipette 50 mL of the prepared sample into a 150 mL plastic beaker. The sample is ready for analysis. The "sample weight" entered in the method parameters is equal to the weight of the sample added to the 500 mL flask multiplied by 0.1. For example, 0.2501 X 0.1 = 0.02501.

**Statistics**

**# of Trials** 5      **Mean** 4962      **%CV** 1.3      **Analysis Time** ~2.5 min

**Comments** 1. The sample's pH has to be in between 4 to 5 before the analysis. Calibrate the electrode with pH buffer 7.00 and 4.01 daily. 2. Analyze the prepared sample within a day (24 hours).  
3. Check the 960 printout for slope (S) and spike recovery. Best results are obtained when the slope is between 25 to 30 and the spike recovery is 98% or better. 4. Rinse the electrodes, stirrer, and dispenser probe between measurements with DI water.

**Method Parameters**

<b>Sample Volume/Weight</b>	0.02501 g	<b>Timed or Stability Readings</b>	3.0 mV/min
<b>Constant Increment</b>	9.0 mV	<b>Number of Endpoints</b>	
<b>Max Titrant Volume</b>	10.0 mL	<b>Desired Units</b>	mg/100g
<b>Molecular weight</b>	40.08	<b>Predose</b>	0.0 mL
<b>Prestir</b>	10.0 sec	<b>Additional Parameters</b>	
<b>Reaction Ratio</b>	1.0		