



**Overview**

The potassium content in potatoes is determined by a multiple known addition technique using 0.098 M potassium standard on the Orion 930 Ionalyzer. A half-cell potassium electrode is used along with a double junction reference electrode, and the 930 Calculates the result and reports it as ppm (w/w).

<b>Market</b>	Food and Beverage	<b>Species Measured</b>	Potassium
<b>Sample</b>	Potato Products	<b>Sample Size</b>	Approximately 10g
		<b>Typical Concentration</b>	~17000 ppm
<b>Technique #</b>	2 Multiple Known Addition	<b>Electrode</b>	Potassium 9319BN; DJ Ref 900200

**Solutions** Potassium Known Addition Standard (Cat. 921906); Potassium Ionic Strength Adjustor (Cat. 931911); Reference Electrode Fill Solution (inner) (Cat. 900002); Reference Electrode Fill Solution (outer) (see electrode manual).  
Thermo Orion 930 (Cat. 093000)

**Sample Prep** For sample D (cooked potatoes) crush and mix several potatoes with a mortar and pestle. Weigh and record approximately 10g of crushed potato, and quantitatively transfer to a blender filled with deionized water to just above the level of the blades. Allow

blender to run at maximum speed for 1 minute, or until is no visible solid left. Transfer solution to a 500 ml volumetric flask using a funnel to insure that no sample is lost. Rinse the cover and inside the blender repeatedly with small aliquots of deion

**Statistics**

**# of Trials** 10      **Mean** 17873 ppm      **%CV** 1.62      **Analysis Time**

**Comments** Rinse the electrodes, stirrer, and dispenser probe thoroughly between measurements with deionized water.  
Potato solution is kept on the magnetic stirrer to make sure samples are homogeneous.  
between samples. Sample weight entered is one tenth of the weighed sample.

**Method Parameters**

<b>Sample Volume/Weight</b>	1.019 g	<b>Timed or Stability Readings</b>	2.0 mv/min
<b>Constant Increment</b>	18.0 mv	<b>Number of Endpoints</b>	
<b>Max Titrant Volume</b>	15.00 ml	<b>Desired Units</b>	ppm-w
<b>Molecular weight</b>	150.2	<b>Predose</b>	0
<b>Prestir</b>		<b>Additional Parameters</b>	Std 0.0980 M; Precision = 2.0 %; TotVol = 51.00 ml
<b>Reaction Ratio</b>	1.0		