

Analytical Performance of Capillary High-Performance Anion-Exchange with Pulsed Amperometric Detection

Jun Cheng, Petr Jandik, Yan Liu, Victor Barreto, and Christopher Pohl
Thermo Fisher Scientific, Sunnyvale, CA, USA

Executive Summary

Capillary-scale high-performance anion-exchange with pulsed amperometric detection (HPAE-PAD) offers substantial benefits over conventional-scale HPAE-PAD, including improved mass detection limits, smaller sample volumes, and reduced eluent consumption. Capillary HPAE-PAD is ideally suited for analytical determination of carbohydrates in a range of matrices across multiple applications.




Figure 1. Block diagram of a capillary Reagent-Free™ HPAE-PAD system with electrochemical detection.

Key Words

Capillary HPAE-PAD, ICS-5000,
Carbohydrates

Abstract

Capillary high-performance anion-exchange with pulsed amperometric detection (HPAE-PAD) is a technique that uses capillary columns (e.g., 0.4 mm i.d. columns) packed with ion-exchange resin and a newly designed capillary amperometric cell with a gold working electrode and a palladium hydrogen (PdH) reference electrode. The cell body is made of titanium and serves as a counter electrode. An electrolytic eluent generator optimized for operation at capillary flow rates is also included in the system.

To characterize the new technique, the authors generated analytical performance data under a number of different conditions based on eluent composition and flow rate. Separation efficiency and concentration detection limits were comparable with those achieved using analytical-scale chromatography. However, the mass detection limits were significantly improved, and the range of linearity of calibration plots was broader than what is typically achieved using analytical-scale chromatography.

Next, the authors used capillary HPAE-PAD for the analysis of carbohydrates in cell cultures, fermentation broths, protein hydrolyzates, beverages, and food samples. Because of its low flow rates (5–10 $\mu\text{L}/\text{min}$), capillary HPAE-PAD offers the convenience of extremely low eluent usage (1 L of water lasts approximately 2 months). An additional advantage is that the system requires only very small injection volumes (0.4 μL).

Experimental

Instrument: Thermo Scientific™ Dionex™ ICS-5000 Ion Chromatography System, including:

- Autosampler with low-volume injector valve (0.40 μ L)
- Dual pump (DP) with flow range from 1 to 100 μ L/min
- Capillary Eluent Generator (CapEG)
- Capillary Thermo Scientific Dionex CR-ATC Continuously Regenerated Anion Trap Column
- Thermo Scientific™ Dionex™ IC Cube™ cartridge with capillary Thermo Scientific™ Dionex™ CarboPac™ PA20 column, 0.40 \times 150 mm
- Capillary ED cell (130 nL dead volume) equipped with PdH reference electrode and disposable gold working electrode (1 mm diameter)

System Control and Data Processing: Thermo Scientific™ Dionex™ Chromeleon™ 6.8 Chromatography Data System software

Chromatographic Conditions: See figures for detail.




Figure 2. Selected molecular structures: sugar alcohols, monosaccharides, and disaccharides.

Table 1. Calibration performance of capillary and analytical systems.

	Linear Range		Correlation Coefficients	
	Capillary	Analytical	Capillary	Analytical
Fuc	0.024–50	0.023–20	0.9958	0.9932
GalN	0.018–25	0.014–20	0.9998	0.9920
GlcN	0.029–25	0.023–20	0.9949	0.9987
Gal	0.054–25	0.029–20	0.9964	0.9932
Glc	0.056–50	0.027–50	0.9986	0.9961
Man	0.068–50	0.059–50	0.9956	0.9975

1. Reference electrode: PdH (capillary); Ag/AgCl (analytical)

2. Eluent: EG generated (capillary); manually prepared eluent (analytical)

3. Injection volume: 0.40 μ L (capillary); 25 μ L (analytical)

4. Analytical system conditions: flow rate 0.25 mL/min; separation column Dionex CarboPac PA20, 3 \times 150 mm; eluent 10 mM NaOH

Table 2. Limits of detection using capillary and analytical systems.

	LOD (μ M)		LOD (picogram)	
	Capillary	Analytical	Capillary	Analytical
Fuc	0.024	0.023	1.6	94.4
GalN	0.018	0.014	1.3	62.7
GlcN	0.029	0.023	2.1	103.0
Gal	0.054	0.029	3.9	130.9
Glc	0.056	0.027	4.0	121.6
Man	0.068	0.059	4.9	265.7

1. Reference electrode: PdH (capillary); Ag/AgCl (analytical)

2. Eluent: EG generated (capillary); manually prepared eluent (analytical)

3. Injection volume: 0.40 μ L (capillary); 25 μ L (analytical)

4. Analytical system conditions: flow rate 0.25 mL/min; separation column Dionex CarboPac PA20, 3 \times 150 mm; eluent 10 mM NaOH




Figure 3. Response stability over two-week period: six monosaccharides, with PdH reference electrode.




Figure 4. Hydrolysis of a monoclonal antibody.




Figure 5. Separation of monoclonal antibody hydrolysate with and without Thermo Scientific™ Dionex™ AminoTrap™ column. The Dionex AminoTrap column is designed to eliminate amino acid interference. It is installed between injector and separation column.




Figure 6. Preparation of a soluble coffee sample for free carbohydrate analysis.




Figure 7. Separation of coffee sugars.




Figure 8. Separation of wood sugars. Two critical pairs—rhamnose/arabinose and xylose/mannose are separated in a single run.




Figure 9. Separation of cell culture medium yeast extract; peptone-dextrose (YPD).




Figure 10. Separation of sugars, amino sugars, and amino acids in wine sample and standard.




Figure 11. Separation of sugars in juice samples and standard.




Figure 12. Separation of sugar alcohols, mono- and disaccharides found in dietary fiber.




Figure 13. Separation of sugars in chocolate syrup.

Conclusion

Capillary-scale HPAE-PAD using the Dionex ICS-5000 system equipped with a capillary ED cell and PdH reference electrode enables determination of carbohydrates in a range of applications. This capillary-scale HPAE-PAD platform provides analytical performance comparable to conventional-scale HPAE-PAD using a regular ED cell and Ag/AgCl reference electrode, with the additional advantages of minimal eluent consumption (<15 mL/day), extremely small sample sizes (400 nL injection volumes), and improved mass detection limits.

www.thermofisher.com/dionex

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Australia +61 3 9757 4486
Austria +43 1 616 51 25
Benelux +31 20 683 9768
 +32 3 353 42 94
Brazil +55 11 3731 5140

China +852 2428 3282
Denmark +45 36 36 90 90
France +33 1 39 30 01 10
Germany +49 6126 991 0
India +91 22 2764 2735

Ireland +353 1 644 0064
Italy +39 02 51 62 1267
Japan +81 6 6885 1213
Korea +82 2 3420 8600
Singapore +65 6289 1190
Sweden +46 8 473 3380

Switzerland +41 62 205 9966
Taiwan +886 2 8751 6655
UK +44 1276 691722
USA and Canada +847 295 7500

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific