

# High Resolution LC-MS Separation and Characterization of Chemoenzymatic Site-specific Engineered Antibody-Drug Conjugates (ADCs)

Shanhua Lin<sup>1</sup>, Terry Zhang<sup>2</sup>, Jonathan Josephs<sup>2</sup>, Brian Agnew<sup>3</sup>, and Xiaodong Liu<sup>1</sup>

1: Thermo Fisher Scientific, Sunnyvale, CA; 2: Thermo Fisher Scientific, San Jose, CA; 3: Thermo Fisher Scientific, Eugene, OR

## ABSTRACT

**Purpose:** Demonstrate application of the Thermo Scientific™ SiteClick™ labeling technology to the production antibody-drug conjugates (ADCs) and a number of different detection probes.

**Methods:** ADCs were made from SiteClick™ enzyme-based N-glycan labeling of bevacizumab and panitumumab with DIBO-MMAE and DIBO-DFO. F(ab)<sub>2</sub> and scFc fragments from intact, azide-labeled intermediate, and ADCs were analyzed using Thermo Scientific™ MAbPac™ RP column coupled with Thermo Scientific™ Orbitrap™ Fusion™ Lumos™ Tribrid™ mass spectrometer.

**Results:** Enzymatic pre-treatment, or “trimming” of the antibody glycans with GlycINATOR® endoglycosidase S2, prior to azide activation, results in highly-reproducible degrees of labeling (DOLs) of 2.0.

## INTRODUCTION

Antibody drug conjugates (ADCs) are becoming a powerful addition to the monoclonal antibody drug therapies. From a safety and efficacy perspective, ensuring the integrity of antibody-cytotoxin conjugate during drug development and production is important. ADCs created using standard lysine and cysteine conjugation methods have multiple Drug-to-Antibody (DAR) forms. Chromatographic separation of multiple DAR forms is challenging due to the sample heterogeneity. Among the several types of conjugation chemistries, enzyme-based site-specific modification shows great potential by eliminating the interruption of antibody-antigen interaction, and providing a highly reproducible and modular conjugation system when compared to standard lysine and cysteine conjugation.

## MATERIALS AND METHODS

### Enzyme Labeling

In a one-pot activation reaction, antibodies were incubated for 1 hours at 37 °C with GlycINATOR (Genovis), then UDP-GalNAz substrate, β-Gal-T1(Y289L) enzyme, and MnCl<sub>2</sub> were added, and the mixture was incubated overnight (16 hrs) in non-phosphate buffer (pH 7.4) at 30 °C. Excess UDP-GalNAz and β-Gal-T1(Y289L) were removed using 50 kD MW cut-off Amicon ULTRA spin filters. The purified azide-activated antibodies were labeled with DIBO-MMAE, or DIBO-DFO at 20 °C (Scheme 1).

### Generation of scFc and F(ab)<sub>2</sub> Fragments

Antibodies were cleaved using FragIT™ MicroSpin columns (Genovis) to generate scFc and F(ab)<sub>2</sub> fragments. All antibody samples, cleaved and non-cleaved, were exhaustively dialyzed into 50 mM ammonium acetate buffer (pH 6.0) and then dried in the SpeedVac for downstream MS analyses.


### LC-MS analysis

Intact and fragment samples were separated and analyzed on Thermo Scientific™ Vanquish™ UHPLC system coupled to Orbitrap Fusion™ Lumos Tribrid mass spectrometer, using a MAbPac™ RP analytical column, 4.0 μm, 2.1 × 50 mm column (p/n 086648) and H<sub>2</sub>O/FA/acetonitrile mobile phase at 80 °C. The MS acquisition method was set with a full scan at both 15,000 (FWHM, at m/z 200) and 120,000 resolution in positive mode. The method parameters were: AGC 2e5, IT 200 ms, in-source CID 0 eV and 35 eV, scan range: 800-3000, 1000-3500 m/z, spray voltage 3.8 kV, sheath gas 60, aux gas 20, capillary temperature 350 °C, s-lens 30, probe heater temperature 150 °C.

### Data processing

Average mass of scFc and F(ab)<sub>2</sub> fragments from unmodified antibody, azide-labeled intermediate, and ADCs were analyzed by Protein Deconvolution 4.0 using Xtract algorithm for isotopically resolved MS spectra and Respect algorithm for isotopically unresolved MS spectra.

### SCHEME 1. SiteClick™ site-specific labeling of antibodies, conjugation to the chitobiose core



## RESULTS

Bevacizumab (IgG1) and panitumumab (IgG2) were used as the model antibodies to synthesize ADCs with DIBO-MMAE and DIBO-DFO. There are 2 N-glycans per antibody Fc domain, one on each heavy chain. After “trimming” of the antibody glycans with GlycINATOR, there is one labeling site available per heavy chain. Therefore, up to 2 MMAE (or DFO) molecules can be conjugated to 1 mAb molecule. HPLC analyses of bevacizumab (Figure 1) and panitumumab (Figure 2) scFc and F(ab)<sub>2</sub> fragments before and after GlycINATOR glycan cleavage and SiteClick labeling shows complete conversion of labeled scFc domains.

HPLC profile of the scFc fragment showed that there are “twin peaks” in each case, regardless of the mAb (bevacizumab or panitumumab) or conjugated molecules (MMAE or DFO). High resolution LC/MS analysis demonstrated that these “twin peaks” have identical masses (Figure 3), indicating that they are structural isomers. NMR studies have shown that DIBO enantiomers were formed after the conjugation.

MS spectra of scFc fragments were acquired using 120K resolution and therefore the deconvolution were carried out using Xtract algorithm which is designed for isotopically resolved peaks. The deconvoluted spectra of scFc and its derivatives (treated with GlycINATOR, labeled with azide, and conjugated to MMAE or DFO) are shown in Figure 4. MS spectra of F(ab)<sub>2</sub> fragments were acquired using 15K resolution and therefore the deconvolution was carried out using ReSpect algorithm which is designed for isotopically unresolved peaks. The deconvoluted spectra of F(ab)<sub>2</sub> are shown in Figure 5. After GlycINATOR cleavage and SiteClick labeling, bevacizumab and panitumumab F(ab)<sub>2</sub> fragments remain completely unaltered ensuring preservation of antigen binding.

Figure 1. HPLC analyses of bevacizumab scFc and F(ab)<sub>2</sub> fragments before and after GlycINATOR glycan cleavage and SiteClick labeling




Figure 2. HPLC analyses of panitumumab scFc and F(ab)<sub>2</sub> fragments before and after GlycINATOR glycan cleavage and SiteClick labeling




Figure 3. LC-MS analyses of bevacizumab-MMAE scFc and F(ab)<sub>2</sub> fragments. The top panel shows the total ion chromatogram (TIC). The MS spectra of scFc fragments eluting at 5.12 min (bottom left) and 5.30 min (bottom right) are shown in the bottom panels.




Figure 4. LC-MS analyses of bevacizumab and panitumumab scFc fragments before and after GlycINATOR cleavage and SiteClick labeling.






Figure 5. LC-MS analyses of panitumumab and panitumumab scFc fragments before and after GlycINATOR cleavage and SiteClick labeling.



## CONCLUSIONS

- SiteClick™, a novel site-specific antibody labeling technology, can be applied to essentially any existing antibody with Fc-domain N-linked glycans.
- This antibody labeling technology has universal application and can be applied to multiple workflows from conjugating antibodies with small organic molecules to large fluorescent probes.
- An additional labeling approach using enzymatic pre-treatment, or “trimming” of the antibody glycans with GlycINATOR, prior to azide activation, enables the chemoenzymatic synthesis of homogeneous antibody conjugates with a degree of labeling of 2 labels per antibody.
- High resolution separation of mAb fragments scFc, scFc-MMAE, scFc-DFO, and F(ab)<sub>2</sub> is achieved on MAbPac RP column within 10 min.
- High resolution MS analysis of scFc fragments and its derivatives leads to successful characterization of site-specific ADCs.

## TRADEMARKS/LICENSES

© 2017 Thermo Fisher Scientific Inc. All rights reserved. GlycINATOR and FragIT are trademarks of Genovis AB. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.