

# Optimization of crosslinked peptide analysis on an Orbitrap Fusion Lumos mass spectrometer

Ryan Bomgarden<sup>1</sup>, Erum Raja<sup>1</sup>, Chris Etienne<sup>1</sup>; Fan Liu<sup>4</sup>, Albert Heck<sup>4</sup>, Mathias Mueller<sup>2</sup>, Rosa Viner<sup>3</sup>

<sup>1</sup>Thermo Fisher Scientific, Rockford, IL; <sup>2</sup>Thermo Fisher Scientific, Bremen, Germany; <sup>3</sup>Thermo Fisher Scientific, San Jose, CA; <sup>4</sup>Utrecht University, the Netherlands

## ABSTRACT

**Purpose:** To improve identification of intra- and inter-protein interactions through analysis of chemically crosslinked peptides.

**Methods:** Different amine-reactive, homobifunctional crosslinkers including disuccinimidyl suberate (DSS), bis-sulfosuccinimidyl suberate (BS3), disuccinimidyl sulfoxide (DSSO)<sup>1</sup> and disuccinimidyl dibutyric urea (i.e. NHS-BuUrBu-NHS)<sup>2</sup> (Figure 1) were compared for protein crosslinking labeling efficiency and crosslinked peptide identification using MS2 and MS3 fragmentation methods. A Thermo Scientific™ Orbitrap™ Fusion Lumos™ mass spectrometer was used for crosslinked peptides analysis. Data analysis was performed by Thermo Scientific™ Proteome Discoverer™ using a XlinkX<sup>3</sup> software node.

**Results:** For both DSSO and BuUrBu, we identified over 40 BSA inter-crosslinked peptides using MS2-MS3 approach compared to less than 20 using MS2 CID for DSSO. We also compared these crosslinkers using an *E. coli* whole cell lysate. Our results show an increase number of identified peptides after crosslinking using the MS2-MS3 in combination with EThcD method compared to CID/EThcD MS2 method.

## INTRODUCTION

Chemical crosslinking in combination with mass spectrometry is a powerful method to determine protein-protein interactions. This method has been applied to recombinant and native protein complexes, and more recently, to whole cell lysates or intact unicellular organisms in efforts to identify protein-protein interactions on a global scale. In this study, we evaluated traditional non-cleavable and MS-cleavable crosslinkers for crosslinked peptide analysis using an Orbitrap Fusion Lumos mass spectrometer. For MS-cleavable crosslinkers, we also compared different types of fragmentation (CID, ETD) and levels of tandem mass spectrometry (MS2 vs MS3). Our data provided insight to the relative performance of different crosslinking compounds and acquisition parameters relevant for improving identification of protein-protein interaction sites.

## MATERIALS AND METHODS

### Sample Preparation

BS3, DSS, DSSO and BuUrBu were used to crosslink 2mg/ml BSA solubilized in 50mM HEPES pH 8 for 1hr at various molar excess of crosslinker to protein. After crosslinking, reactions were quenched with 1M Tris pH 8 and analyzed by SDS-PAGE or reduced, alkylated and digested with trypsin for MS analysis. Protein and peptide concentrations were determined using the Pierce™ BCA Protein Assay Kit and the Pierce™ Quantitative Colorimetric Peptide Assay, respectively. *E. coli* lysates were crosslinked using 20-fold molar excess of crosslinker to protein before reduction, alkylation and digestion. Peptides were fractionated using a HyperSep™ Retain CX column (30mg) with an increasing step gradient of ammonium acetate (e.g. 50mM, 150mM, 250mM, 250mM, 500mM, 1M). Fractionated samples were desalted using C18 before LC-MS/MS analysis.

### Liquid Chromatography and Mass Spectrometry

Samples were separated by RP-HPLC using a Dionex™ Ultimate 3000 system connected to a Thermo Scientific™ EASY-Spray™ column, 50 cm × 75 μm over a 1 hr. 4-40% gradient (A: water, 0.1% formic acid; B: acetonitrile, 0.1% formic acid) at 300 nL/min flow rate. The crosslinked BSA and *E. coli* cell lysates samples were analyzed on Thermo Scientific™ Orbitrap Fusion™ Lumos™ and QExactive™ HF mass spectrometers. Additional LC and MS settings are shown in Table 1.

### Data Analysis

Spectral data files were analyzed by XlinkX 2.0 or Thermo Scientific™ Proteome Discoverer™ 2.2 software using the XlinkX node for crosslinked peptides and SEQUEST@HT search engine for unmodified and dead-end-modified peptides. Carbamidomethylation (+57.021 Da) used as a static modification for cysteine. Different crosslinked mass modifications for lysine were used as variable modifications for lysine in addition to methionine oxidation (+15.996 Da). Data was searched against a Swiss-Prot® *E. coli* or BSA databases with a 1% FDR criteria for protein spectral matches. For MS2-MS3 methods, a linear-peptide search option (using MS3 scans for identification and MS2 scan for detection of crosslinked peptides) was used for XlinkX database searching. The XlinkX standard enumeration search option was used for data acquired using the MS2 methods (e.g. CID, ETD, EThcD).<sup>3</sup>

Figure 1. Structures of non-cleavable and MS-cleavable crosslinkers used for protein-protein interaction analysis.





Table 1. LC-MS acquisition and database search parameter settings.

	Orbitrap Fusion Lumos Parameters	Q Exactive HF Parameters
<b>LC gradient</b>	6-40% in 45min	4-40% in 65min
<b>Full MS</b>	OT	OT
<b>Resolution</b>	120K	120K
<b>Target value</b>	2e5	3e6
<b>Max injection time</b>	100	50
<b>Top N</b>	5 sec	15
<b>MS2</b>	OT CID	OT HCD
<b>Isolation mode</b>	Quadrupole	Quadrupole
<b>Isolation width</b>	1.6	1.4
<b>NCE</b>	25	30, SID 15-25
<b>Resolution</b>	30K	15K
<b>Target value</b>	5e4	1e5
<b>Max injection time</b>	100ms	100ms
<b>MS3 SPS</b>	IT HCD	
<b>Isolation width</b>	2	
<b>NCE</b>	30	
<b>Resolution</b>	Rapid	
<b>Target value</b>	2e4	
<b>Max injection time</b>	120 ms	
<b>Search parameters</b>	XlinkX, SequestHT	XlinkX, SequestHT
<b>Precursor tolerance</b>	10ppm	10ppm
<b>Fragment tolerance</b>	0.02Da	0.02Da
<b>Static</b>	Carbamidomethyl (C)	Carbamidomethyl (C)
<b>Dynamic</b>	Oxidation (M) DSS, DSSO or BuUrBu(K)	Oxidation (M) DSS, DSSO or BuUrBu (K)

Figure 2. MS acquisition used for MS2 (A) or MS2-MS3 (B) fragmentation methods.



## RESULTS

Figure 3 Comparison of BSA crosslinking efficiency by SDS-PAGE. Different crosslinkers were incubated with BSA at molar excess of crosslinker to protein (e.g. 20, 100 or 500-fold). Crosslinking efficiency is shown by decreased mobility by SDS-PAGE and varied by crosslinker type, solubility and concentration.




Figure 4 BSA crosslinked peptide spectra identified by MS2-MS3 method and XlinkX using DSSO (A) or BuUrBu (B) crosslinkers. XlinkX uses unique fragment ion patterns of MS-cleavable crosslinkers (purple annotation) to detect and filter crosslinked peptides for a crosslinked database search.




Figure 5. Graph showing number of BSA crosslinked peptides identified using different non-cleavable (BS3, DSS) and cleavable crosslinkers (DSSO, BuUrBu) for various MS<sup>n</sup> methods. Both BS3 and DSS had similar numbers of crosslinked peptides identified for CID and EThcD methods. BuUrBu had more crosslinked peptides identified by CID and SID-HCD on a QExactive. Although DSSO, had the fewest crosslinked peptides identified by CID and HCD, it had the most for the MS2-MS3 method if the linear-peptide search mode<sup>3</sup> is used. All crosslinkers showed similar numbers of identified crosslinked peptides by EThcD.




Figure 6. The processing (A) and consensus (B) XlinkX workflows in Proteome Discoverer 2.2 software including a separate crosslinkers results tab (C) and spectra annotation (D).




Figure 7. Comparison of different EThcD energies for crosslinked peptide fragmentation. Increasing EThcD fragmentation energy results in different fragment ion intensity in MS/MS spectra (A) and unique identified crosslinked peptides (B).





Figure 8. *E. coli* cell lysate crosslinked peptides identified using different instrument methods and XlinkX software. MS2-MS3 method combined with MS2 EThcD (A) provided most identifications (B) compared to MS2 CID, MS2 CID-MS3 ETD, or MS2 CID-MS3 HCD methods.



## CONCLUSIONS

- MS-cleavable crosslinkers, DSSO and BuUrBu, crosslink BSA with slightly lower efficiency than non-cleavable crosslinkers, DSS and BS3, possibly due to small differences in crosslinker length or solubility.
- DSS, BS3 and BuUrBu worked well for CID, HCD and EThcD MS2 fragmentation methods. However, DSSO resulted in the most identified BSA crosslinked peptides using the a combination of MS2, MS3 spectral sequence information and XlinkX.
- Different EThcD energies not only changed MS/MS fragment ion intensities but also resulted in different identified crosslinked peptides.
- For more complex *E. coli* crosslinked samples, using the MS2-MS3 in combination with EThcD method resulted the most identified crosslinked peptides compared to other methods.

## REFERENCES

- Kao A, Chiu CL, Vellucci D, Yang Y, Patel VR, Guan S, Randall A, Baldi P, Rychnovsky SD, Huang L. Development of a novel cross-linking strategy for fast and accurate identification of cross-linked peptides of protein complexes. 2011. *Mol Cell Proteomics*, 10(1): M110.002212.
- Müller MQ, Dreierock F, Ihling CH, Schäfer M, Sinz A. Cleavable cross-linker for protein structure analysis: reliable identification of cross-linking products by tandem MS. 2010. *Anal Chem*. 82(16):6958-68.
- Liu F, Rijkers D, Post H, Heck AJ. Proteome-wide profiling of protein assemblies by cross-linking mass spectrometry. 2015. *Nat Methods*, 12(12):1179-84.

## ACKNOWLEDGEMENTS

The authors would like to thank Kai Fritzscheier (Thermo Fisher Scientific, Germany) and Richard Schellema (University of Utrecht) for their work on integrating XlinkX software as a node in Proteome Discoverer 2.2.

## TRADEMARKS/LICENSING

For Research Use Only. Not for use in diagnostic procedures.  
© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.