

# Simultaneous N, C and S Isotope Ratio Determination on a DELTA V Isotope Ratio MS using a Flash Elemental Analyzer

Oliver Kracht, Thermo Fisher Scientific, Bremen, Germany

## Introduction

Multi-element isotope ratio analysis is becoming increasingly important, particularly in food authenticity control and forensic applications. With increasing numbers of samples for analysis, reduced sample preparation and measurement time are required. The ability to measure N, C and S isotopes from a single sample drop was demonstrated by Avak in 1996 with an analysis time of 30 minutes.<sup>1</sup>

With the introduction of the Thermo Scientific ConFlo IV universal interface in July 2007, the ability for multi-element isotope ratio determination was both enhanced and simplified. The ConFlo IV interface allows software-controlled switching of up to five reference gases and automatic dilution of sample gases down to 1:100 according to the varying amounts of the three elements in the samples. Using the Thermo Scientific smartEA option, the simultaneous analysis of N, C and S is simplified even further. The smartEA™ option adjusts the dilution of sample gases as a direct response to the TCD signals from the Flash Elemental Analyzer (EA). The optimized chromatography now reduces the analysis time to less than 18 minutes.

This technical note explains the modifications needed for triple measurements of N, C and S isotope ratios and shows the performance that can be achieved.

## Method

Optimal conditions for sulfur isotope ratio analysis require a short transfer line and reduced volumes in a water-free environment. Therefore a single reactor filled with WO<sub>3</sub> and electrolytic copper wires is used for the triple N, C and S analysis. A chemical trap of magnesium perchlorate removes water produced in the combustion process. Tubing in the Flash EA and the connection to the ConFlo IV interface is made of 1/16" Sulfinert® capillary to avoid the formation of acidic compounds by SO<sub>2</sub> with adsorbed water. The commonly used Teflon® tubing cannot be used for triple analysis of N, C and S due to increased nitrogen background. An optimized chromatographic column separates N<sub>2</sub> and CO<sub>2</sub> with an improved retention time of the SO<sub>2</sub>. The Flash elemental analyzer used was a Thermo Scientific FlashEA 1112, system parameters are given in Table 1.




Figure 1: N, C and S triple measurement of sulfanilamide. Top: mass traces, middle: TCD trace of Flash EA, bottom: tabs with comprehensive data and results for each isotope.

### FlashEA 1112 settings NCS triple analysis

|                     |           |
|---------------------|-----------|
| Reactor temperature | 1,020 °C  |
| GC temperature      | 70 °C     |
| Carrier flow        | 80 mL/min |
| Autosampler type    | MAS 200R  |

Table 1: System settings for triple measurement.

The ConFlo IV universal interface and the Flash EA were connected with the smartEA™ option to make use of the automatic dilution ability. The isotope ratio mass spectrometer (IRMS) used was a Thermo Scientific DELTA V Plus with a universal triple collector. A fast magnet field jump of the IRMS takes place automatically between the N<sub>2</sub> and CO<sub>2</sub> peaks followed by a slow magnet field jump to SO<sub>2</sub>.

## Results

The chromatogram in Figure 1 shows flat top reference gas pulses from the ConFlo IV interface and the chromatographically separated sample peaks of N<sub>2</sub>, CO<sub>2</sub> and SO<sub>2</sub> from the elemental analyzer (400 µg of sulfanilamide). N<sub>2</sub> and CO<sub>2</sub> show a separation of 17 seconds which is required to perform a fast magnet field jump immediately after the detection of the N<sub>2</sub> sample peak at 191 seconds. The change of the magnet field from CO<sub>2</sub> to SO<sub>2</sub> is performed slowly between the CO<sub>2</sub> and SO<sub>2</sub> reference gas pulses at 460 seconds.

## Key Words

- ConFlo IV
- DELTA V Series
- Elemental Analyzer
- smartEA™
- Multi-Element Analysis
- Isotope Ratio MS

| Sample Identifier | Amount [mg] | $\delta^{15}\text{N}$ [‰] | $\delta^{13}\text{C}$ [‰] | $\delta^{34}\text{S}$ [‰] | Mean $\delta^{15}\text{N}$ [‰]<br>S.D. | Mean $\delta^{13}\text{C}$ [‰]<br>S.D. | Mean $\delta^{34}\text{S}$ [‰]<br>S.D. | C/S Ratio |
|-------------------|-------------|---------------------------|---------------------------|---------------------------|--|--|--|-----------|
| Sulfanilamide     | 0.245       | -1.76                     | -28.76                    | 3.09                      | -1.84                                  | -28.77                                 | 3.04                                   | 2.2       |
| Sulfanilamide     | 0.192       | -1.87                     | -28.80                    | 2.90                      | <b>0.07</b>                            | <b>0.03</b>                            | <b>0.12</b>                            |           |
| Sulfanilamide     | 0.244       | -1.89                     | -28.73                    | 3.13                      |  |  |  |           |
| Vitamin B1        | 0.192       | -1.02                     | -32.88                    | -4.94                     |  |  |  | -0.96     |
| Vitamin B1        | 0.240       | -0.89                     | -32.86                    | -5.06                     | <b>0.06</b>                            | <b>0.04</b>                            | <b>0.09</b>                            |           |
| Vitamin B1        | 0.196       | -0.96                     | -32.81                    | -4.89                     |  |  |  |           |
| Vitamin H         | 0.412       | -12.34                    | -29.23                    | 10.42                     |  |  |  | -12.30    |
| Vitamin H         | 0.468       | -12.31                    | -29.30                    | 10.32                     | <b>0.04</b>                            | <b>0.03</b>                            | <b>0.10</b>                            |           |
| Vitamin H         | 0.474       | -12.27                    | -29.24                    | 10.52                     |  |  |  |           |
| Cysteine          | 0.299       | 10.03                     | -14.26                    | 7.13                      |  |  |  | 10.12     |
| Cysteine          | 0.288       | 10.32                     | -14.23                    | 7.17                      | <b>0.17</b>                            | <b>0.02</b>                            | <b>0.07</b>                            |           |
| Cysteine          | 0.280       | 10.01                     | -14.27                    | 7.03                      |  |  |  |           |
| Soil              | 95.540      | 7.01                      | -23.61                    | 3.73                      |  |  |  | 7.00      |
| Soil              | 100.979     | 7.05                      | -23.67                    | 4.23                      | <b>0.05</b>                            | <b>0.05</b>                            | <b>0.25</b>                            |           |
| Soil              | 103.211     | 6.95                      | -23.57                    | 4.04                      |  |  |  |           |
| Octopus Tissue    | 2.041       | 10.48                     | -16.09                    | 15.91                     |  |  |  | 10.50     |
| Octopus Tissue    | 2.033       | 10.49                     | -16.06                    | 16.10                     | <b>0.02</b>                            | <b>0.02</b>                            | <b>0.11</b>                            |           |
| Octopus Tissue    | 2.070       | 10.52                     | -16.06                    | 15.91                     |  |  |  |           |
| Bird Feather      | 1.148       | 7.60                      | -15.42                    | 5.42                      |  |  |  | 7.71      |
| Bird Feather      | 1.002       | 7.75                      | -15.58                    | 5.13                      | <b>0.10</b>                            | <b>0.09</b>                            | <b>0.17</b>                            |           |
| Bird Feather      | 1.065       | 7.78                      | -15.57                    | 5.14                      |  |  |  |           |
| Pesticide         | 0.445       | -9.28                     | -31.78                    | -4.77                     |  |  |  | -9.21     |
| Pesticide         | 0.457       | -9.01                     | -31.77                    | -4.73                     | <b>0.17</b>                            | <b>0.06</b>                            | <b>0.11</b>                            |           |
| Pesticide         | 0.438       | -9.34                     | -31.67                    | -4.56                     |  |  |  |           |
| River Sediment    | 5.173       | 6.00                      | -22.75                    | -22.65                    |  |  |  | 6.37      |
| River Sediment    | 5.327       | 6.65                      | -22.42                    | -21.89                    | <b>0.33</b>                            | <b>0.17</b>                            | <b>0.38</b>                            |           |
| River Sediment    | 5.304       | 6.46                      | -22.52                    | -22.36                    |  |  |  |           |

Table 2: Various sample types analyzed simultaneously for N, C and S isotope ratios. Samples were measured in triplicates.

The smartEA option delays the peaks after TCD detection. The shift in retention time of the sample peaks between TCD and mass trace is 60 seconds. This time is needed to adjust the automatic sample dilutions in the ConFlo IV interface. Each gas peak is individually diluted with respect to its abundance in the sample. This results in an acquisition time of 19 minutes which allows the analysis of three samples per hour.

The precision (S.D.) shown in Table 2 is  $\leq 0.17\%$  for  $\delta^{15}\text{N}$ ,  $\leq 0.1\%$  for  $\delta^{13}\text{C}$  and  $\leq 0.17\%$  for  $\delta^{34}\text{S}$ . For heterogeneous samples and samples with low sulfur contents (e.g. soil, river sediment, see Table 2) reduced precision has to be expected. The C/S ratio should not exceed 70. Sulfur contents are best analyzed from 30  $\mu\text{g}$  on.

## Conclusion

The triple analysis of N, C and S isotopes can be performed with any DELTA V series mass spectrometer using a Flash elemental analyzer for IRMS with a sulfur applicable reactor, tubing, and a chromatographic column described in Table 3. It is recommended to use Sulfinert capillaries instead of permeable Teflon or stainless steel tubing in order to reduce background and to avoid water adsorption. The analysis time can be reduced to less than 18 minutes.

| Parts   | Possible Supplier        | Part Number |
|---|--------------------------|-------------|
| NCS Reactor (sulfur reactor with $\text{WO}_3$ and reduced Cu)    | Thermo Fisher Scientific | 468 020 21  |
| Adsorption Trap (Water Trap) with connectors                      | Thermo Fisher Scientific | 281 131 02  |
| Aluminium Ferrule, 1.6 mm (Olive, 10 pcs)                         | Thermo Fisher Scientific | 290 340 44  |
| Stainless Steel Tubing Nut ss 304 6MB (10 pcs)                    | Thermo Fisher Scientific | 350 404 03  |
| Sulfinert Separation Column PP-QS, 2 m, 1/4", 50/80 mesh          | Restek                   | PC4974      |
| Sulfinert Tubing (1.5 m plus distance EA-ConFlo) 0.03"ID, 1/16"OD | Restek                   | 22504       |

Table 3: Required parts for simultaneous N, C and S triple isotope analysis.

This setup can also be used to analyze N or C or S in a single mode. Analysis time for single sulfur isotope analysis can be reduced by increasing the column temperature allowing for a faster  $\text{SO}_2$  elution to the disadvantage of the  $\text{N}_2/\text{CO}_2$  separation.

Modifications and maintenance are minor and equal to a dedicated setup for sulfur isotope analysis. All parts can be installed quickly and easily on site at reasonable costs.

## References

- Avak, H. Triple Element Mode – C, N, S from a single combustion. *Application Flash Report No. 16*, Finnigan MAT, 1996.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

**Africa-Other**  
+27 11 570 1840

**Australia**  
+61 3 9757 4300

**Austria**  
+43 1 333 50 34 0

**Belgium**  
+32 53 73 42 41

**Canada**  
+1 800 530 8447

**China**  
+86 10 8419 3588

**Denmark**  
+45 70 23 62 60

**Europe-Other**  
+43 1 333 50 34 0

**Finland/Norway/Sweden**  
+46 8 556 468 00

**France**  
+33 1 60 92 48 00

**Germany**  
+49 6103 408 1014

**India**  
+91 22 6742 9434

**Italy**  
+39 02 950 591

**Japan**  
+81 45 453 9100

**Latin America**  
+1 561 688 8700

**Middle East**  
+43 1 333 50 34 0

**Netherlands**  
+31 76 579 55 55

**New Zealand**  
+64 9 980 6700

**Russia/CIS**  
+43 1 333 50 34 0

**South Africa**  
+27 11 570 1840

**Spain**  
+34 914 845 965

**Switzerland**  
+41 61 716 77 00

**UK**  
+44 1442 233555

**USA**  
+1 800 532 4752

[www.thermoscientific.com/irms](http://www.thermoscientific.com/irms)

©2010, 2011 Thermo Fisher Scientific Inc. All rights reserved. Sulfinert is a registered trademark of Restek Inc. Teflon is a registered trademark of Dupont Inc. ISO is the registered trademark of the International Standards Organisation All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

AN30194\_E 11/11G