

Gain Calibration Protocol for $10^{13} \Omega$ Resistor Current Amplifiers using the Certified Neodymium Standard JNdi-1 on the TRITON *Plus*

Anne Trinquier, Thermo Fisher Scientific, Bremen, Germany

Key Words

$10^{13} \Omega$ current amplifiers, $10^{11} \Omega$ current amplifiers, gain calibration, TRITON *Plus*, thermal ionization MS, static mode, neodymium, JNdi-1

Introduction

Newly developed $10^{13} \Omega$ current amplifiers,¹ allow analysis of minor ion beam intensities in the range 3aA ($3 \cdot 10^{-18}$ A, equivalent to 2000 cps on a ion counter) to 0.5 fA ($5 \cdot 10^{-16}$ A) on Faraday cups. They can be used alone when all ion beam intensities are low, or in combination with $10^{11} \Omega$, $10^{12} \Omega$ or $10^{10} \Omega$ current amplifiers, for measuring the ion beams of the minor isotopes from elements with extreme isotopic ratios (e.g. Ca, U, Sr) and/or for simultaneously collecting isobaric interfering species, thus providing precise and accurate isobaric interference correction.

Selective virtual amplifier concept²

Following the virtual amplifier concept, during analysis, all cups are sequentially connected to all active amplifiers, all ion beams are in turn measured with the same set of amplifiers, and the stochastic calibration biases of the amplifiers can be averaged out. The virtual amplifier thus reduces the propagation of the uncertainty of the gain calibration procedure. With the introduction of high resistor current amplifiers (10^{12} and $10^{13} \Omega$ current amplifiers), a selective virtual amplifier tool has been developed that allows sequentially connecting all active $10^{11} \Omega$ amplifiers to all corresponding cups, while maintaining the connection between the 10^{12} and $10^{13} \Omega$ amplifiers to their respective cups. This selective virtual amplifier protocol allows simultaneously high-precision and accurate isotopic ratios assessment (if one takes the example of neodymium, all Nd isotopes⁴⁻⁵) and precise and accurate analysis of minor isotopes of the element of interest and/or isobaric interfering species (in this example: Ce, Sm).

The present technical note presents a gain calibration protocol for $10^{13} \Omega$ resistor current amplifiers using the certified neodymium standard JNdi-1 on the Thermo Scientific™ TRITON *Plus*™ Thermal Ionization MS (TIMS).

Choice of a Reference: Certified Neodymium Standard JNdi-1

Neodymium isotopic ratios obtained for JNdi-1 (Split 1 position 179, ref. 3) over 3 sessions of analyses serve as a reference for this study.³ Neodymium has been chosen because up to 7 Nd isotopes can be measured simultaneously and to comparable ion beam intensities. JNdi-1 has been selected because it is a readily accessible standard with low levels of interfering contaminants (Ce, Sm). Mass bias corrected isotopic ratios $^{142}\text{Nd}/^{144}\text{Nd}$, $^{143}\text{Nd}/^{144}\text{Nd}$, $^{145}\text{Nd}/^{144}\text{Nd}$, $^{148}\text{Nd}/^{144}\text{Nd}$ and $^{150}\text{Nd}/^{144}\text{Nd}$ of JNdi-1 have been measured on $10^{11} \Omega$ current amplifiers using the virtual amplifier, following the protocol previously described⁴⁻⁵ and provide the reference isotopic ratios for this study. Three sessions of analysis⁶ providing 2–5 ppm/amu external reproducibility (2RSD, n=15, 500 ng loads, 9.9 V $^{142}\text{Nd}^+$, July 2013), 1–2 ppm/amu (2RSD, n=3, 300-500 ng loads, 7.4 V $^{142}\text{Nd}^+$, January 2014) and 2–8 ppm/amu (2RSD, n=7, 200 ng loads, 4.6V $^{142}\text{Nd}^+$, February 2014) have been combined. The resulting external reproducibility of Nd isotopic ratios over the 3 sessions of analysis is 2–10 ppm/amu (2RSD, n=25, 200-500 ng loads).

ng	^{142}Nd (V)	$^{142}\text{Nd}/^{144}\text{Nd}$	2RSD	2RSD (ppm)	$^{143}\text{Nd}/^{144}\text{Nd}$	2RSD	2RSD (ppm)	$^{145}\text{Nd}/^{144}\text{Nd}$	2RSD	2RSD (ppm)	$^{148}\text{Nd}/^{144}\text{Nd}$	2RSD	2RSD (ppm)	$^{150}\text{Nd}/^{144}\text{Nd}$	2RSD	2RSD (ppm)
400	8	1.141836	9E-06	8	0.512107	5E-06	9	0.348404	3E-06	10	0.241580	2E-06	9	0.236452	4E-06	17

Figure 1. Neodymium isotopic composition of JNdi-1.

The ratios $^{142}\text{Nd}/^{144}\text{Nd}$ and $^{143}\text{Nd}/^{144}\text{Nd}$ agree within respective uncertainty with literature: $^{142}\text{Nd}/^{144}\text{Nd} = 1.1418351 \pm 42$ (ref.7); $^{143}\text{Nd}/^{144}\text{Nd} = 0.512115 \pm 7$ (ref.3); $^{145}\text{Nd}/^{144}\text{Nd} = 0.512108 \pm 6$ (ref.8). All ratios agree within respective uncertainty with the data from ref. 9.

Line No.	Mass Set	L4	L3	L2	L1	RPQ/IC1 C	H1	H2	H3	H4	Integration Time(s)	Number of Integrations	Idle Time(s)	Control Cup Peakcenter	Control Cup Focus
1	Main	^{140}Ce	^{142}Nd	^{143}Nd	^{144}Nd	^{145}Nd	^{148}Nd	^{147}Sm	^{149}Nd	^{150}Nd	16.77	1	18.000	^{144}Nd	^{144}Nd

Figure 2. Neodymium Cup configuration.

Gain Calibration Protocol for $10^{13} \Omega$ Current Amplifiers

In the present study, $10^{11} \Omega$ current amplifiers are connected to the cups collecting ^{144}Nd and ^{146}Nd . ^{144}Nd is the denominator of all Nd isotopic ratios and $^{144}\text{Nd}/^{146}\text{Nd}$ is the reference ratio used for instrumental mass bias correction, using the exponential law. Up to five $10^{13} \Omega$ current amplifiers can be connected to masses ^{142}Nd , ^{143}Nd , ^{145}Nd , ^{148}Nd , ^{150}Nd , allowing up to five $10^{13} \Omega$ current amplifiers gains to be assessed simultaneously. JNdi-1 loads of 200–500 ng are run at ion beam intensities <500 mV on $^{142}\text{Nd}^+$ for 80 cycles, 17 s integration times and an analysis time of ca 40 min. The ratio between certified and mass bias corrected $^{142}\text{Nd}/^{144}\text{Nd}$ values multiplied by 0.01 (to account for the 2 order of magnitude difference between the 10^{11} and the $10^{13} \Omega$ amplifiers) yields the amplifier gains.

Note: Ion beams measurable with $10^{13} \Omega$ current amplifiers are highlighted in red, ion beams measured with $10^{11} \Omega$ current amplifiers are highlighted in blue. In this study, up to four $10^{13} \Omega$ current amplifiers gains have been assessed simultaneously using masses ^{142}Nd (cup L3), ^{143}Nd (cup L2), ^{145}Nd (cup C), ^{148}Nd (cup H3).

In this campaign of analyses, the gains of two $10^{13} \Omega$ current amplifiers were assessed from the period September 2013 to February 2014 and the gains of two additional $10^{13} \Omega$ current amplifiers were assessed from the period September to November 2013. Internal precision is on average better than 50 ppm (2RSE) and external reproducibility over 6 months is better than 100 ppm (2RSD).




Figure 3. Gain calibration external reproducibility (2RSD) over 6 months (#1-2) and over 3 months (#3-4).

References

1. Bouman *et al.*, Thermo Fisher Scientific Technical Note 30282.
2. Tuttas and Schwieters, Thermo Fisher Scientific Application Note 30015.
3. Tanaka *et al.* (2000) *Chem. Geol.* **168**, 279-281.
4. Trinquier, Thermo Fisher Scientific Application Note 30280.
5. Trinquier, Thermo Fisher Scientific Application Note 30281.
6. Trinquier, Thermo Fisher Scientific Application Note in prep.
7. Roth *et al.* (2013) *EPSL* **361**, 50-57.
8. O'Neil *et al.* (2008) *Science* **321**, 1828-1831.
9. Brennecka *et al.* (2013) *PNAS* **110**, 17241-17246.

thermoscientific.com/irms

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
 400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00
Russia/CIS +43 1 333 50 34 0

Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
 SCIENTIFIC

A Thermo Fisher Scientific Brand