

Separation and Identification of Polar Metabolites by Ion Chromatography Coupled with High-Resolution Accurate-Mass Mass Spectrometry

Monika Verma and Terri Christison
Thermo Fisher Scientific, Sunnyvale, CA, USA

Executive Summary

Analyzing small metabolites is necessary to identify biomarkers of certain disease states. The separation and identification of different classes of metabolites requires multiple analytical techniques. Ion chromatography (IC) is an effective technique for separating polar metabolites as it uses both a polar stationary phase and a polar mobile phase. Additionally, High-Resolution Accurate-Mass (HR/AM) mass spectrometry can confirm the identity of the given metabolite. When IC is coupled with HR/AM mass spectrometry, a large quantity of polar metabolites can successfully be determined at low concentrations. This further complements other analytical techniques used for metabolite analysis including Ultra High Performance Liquid Chromatography (UHPLC), Hydrophilic Interaction Liquid Chromatography (HILIC), and HR/AM mass spectrometry.

Keywords

Metabolomics, isobaric compounds, Ion Chromatography, capillary RFIC, mass spectrometry

Metabolomics Analysis

The study of metabolomic phenotypes quantifies numerous small molecule components called metabolites found in biological fluids and cell extracts. The relative quantitation of an entire complement of low molecular weight metabolites in plant or animal samples is compared during the process known as metabolic profiling. Differences in the up or down regulation of metabolites in a disease state relative to that of a normal state allow for the identification of many potential biomarkers in a given sample. This facilitates identifying potential causes and cures of the disease. To fully understand these metabolomic changes in the disease state, samples are prepared in controlled experiments: controls (non-disease state or non-aggressive disease state), disease state (or aggressive disease state), and samples with gene-modifications associated with the disease (knock-down). These samples typically require lengthy generation and sample preparation times before analysis to separate and identify the metabolomic differences of the control, disease state, and knock-down samples.

The preparation and analysis of these typically small volume samples can be quite challenging. Careful handling and storage of these precious biological samples are important as a large amount of information must be obtained from limited samples. The quality of samples is affected by factors including but not limited to collection time, containers used, and the types of preservatives and additives.

These biological samples also tend to include complex mixtures of cell materials and small organic molecules including but not limited to organic acids, inorganic anions and cations, sugars, and sugar phosphates. In addition, these molecules are found in disparate concentrations within the samples and exhibit a broad range of polarities. The variety of these small metabolites requires the use of more than one analytical technique to effectively separate and identify isobaric compounds.¹

Separation and Identification of Metabolites

Common analytical techniques utilized to separate and identify metabolites includes the use of traditional reversed-phase chromatography, UHPLC or HILIC. Traditional reversed-phase columns using a non-polar stationary phase are most effective for hydrophobic metabolites. The use of UHPLC or HILIC is generally a viable option, but can be limiting for resolving polar compounds. Additionally, many of these compounds are isobaric which can prevent sole identification by mass spectrometry. However, ion chromatography is an alternative technique that provides improved quantification and sensitivity of polar compounds. An application example will be highlighted and discussed later in Figure 3. A typical IC system is shown in Figure 1.




Figure 1. Ion chromatography system.

IC Coupled with High-Resolution Accurate-Mass Mass Spectrometry

When utilized in combination, IC coupled with HR/AM mass spectrometry can successfully separate and identify a larger quantity of polar metabolites at low concentrations compared to UHPLC or HILIC. Direct coupling of IC with HR/AM mass spectrometry is possible with Reagent-Free™ IC (RFIC™) systems as the hydroxide mobile phase (eluent) is converted to pure water by the electrolytic eluent suppressor. RFIC systems have automated eluent generation to electrolytically create the required eluents used for IC applications, thereby increasing consistency and reproducibility.

IC with HR/AM Mass Spectrometry Interface

Continuous in-line desalting

In traditional mass spectrometry, the use of a hydroxide mobile phase (eluent) may be of concern due to its non-volatile and corrosive nature. However, in anion-exchange chromatography, electrolytic eluent suppression converts hydroxide mobile phase back to pure water. This occurs via an ion-exchange substitution of K^+ with H^+ ions. The analyte is simply converted to its acid form.

Figure 2 illustrates the chemistry used in an anion-electrolytic suppressor. The suppressor removes potassium from the eluent and replaces it with hydronium ions formed by the electrolysis of the water. These hydronium ions combine with the hydroxyl ions from the eluent to form water.




Figure 2. Dionex electrolytically regenerated anion suppressor.

Highlighted Application

The introduction of capillary format RFIC systems further improves direct coupling with HR/AM mass spectrometry, especially in metabolomic research applications where sample size is commonly limited. A common approach using capillary RFIC systems is to perform a large-volume direct injection, which is suitable for samples with low levels of matrix ions. For example, a 10 μL injection onto a 0.4 mm i.d. column in a capillary IC system is equivalent to a 1000 μL injection onto a 4 mm i.d. column. Therefore, capillary RFIC systems are beneficial when sample volumes are limited.²

Coupling a capillary RFIC system with a HR/AM mass spectrometry has been demonstrated in profiling metabolic biomarkers (sugar monophosphates) of oral squamous cell carcinoma (OSCC) metastasis in cell lysates.¹ Although sugar phosphates are generally difficult to separate by traditional LC methods, sugar phosphates are routinely separated by IC. Figure 3 demonstrates how a capillary RFIC system coupled to HR/AM mass spectrometry provides the superior resolution and analytical sensitivity of these isobaric, polar metabolites as compared to UHPLC and HILIC. Please note that in the capillary RFIC system to mass spectrometer interface, a polar desolvation solvent (methanol) was used for maximizing sensitivity.




Figure 3. Capillary IC-HR/AM results had higher resolution and higher peak intensities over UHPLC- and HILIC-HR/AM spectroscopy using the Thermo Scientific™ Q Exactive™ Mass Spectrometer. Eleven sugar monophosphates ($m/z - 1 = 259.0225$) were resolved by capillary IC-HR/AM with ~100-fold higher intensities as compared to three sugar phosphates by HILIC-HR/AM.

Conclusion

IC with HR/AM mass spectrometry provides additional sensitivity and resolving abilities needed to identify and quantify polar compounds at lower concentrations in metabolomics analysis. Inline desalting by electrolytic suppressors converts the mass spectrometer-incompatible hydroxide mobile phase (eluent) to pure water, thereby making the flow from the IC compatible with the HR/AM mass spectrometer. Additionally, capillary IC-HR/AM mass spectrometry minimizes the consumption of precious biological samples.

References

1. Wang, J., Christison, T., Misuno, K., Lopez, L., Huhmer, A. F., Huang, Y., Hu, S. Metabolomics Profiling of Anionic Metabolites in Head and Neck Cancer Cells by Capillary Ion Chromatography with Orbitrap Mass Spectrometry. Analytical Chemistry [Online] pubs.acs.org (accessed June 3, 2014).
2. Thermo Fisher Scientific. Thermo Scientific Dionex Reagent-Free Ion Chromatography (RFIC) System Capabilities Brochure. Thermo Scientific P/N BR70667_E, Sunnyvale, CA, 2013.
3. Thermo Fisher Scientific. White Paper 70924: Highlighted Capillary IC Applications. Thermo Scientific P/N WP70924, Sunnyvale, CA, 2013.

www.thermoscientific.com/chromatography

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Fisher Scientific, Sunnyvale, CA
USA is ISO 9001:2008 Certified.

Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Brazil +55 11 3731 5140
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
400 650 5118

Denmark +45 70 23 62 60
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 6 6885 1213
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Singapore +65 6289 1190
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
Taiwan +886 2 8751 6655
UK/Ireland +44 1442 233555
USA +1 800 532 4752

WP71102-EN 0814S

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand