

Reacti-Therm

Sample Derivatization System

This unique product combines heating, stirring and evaporation with unmatched convenience and versatility. The system is comprised of a Thermo Scientific™ Reacti-Therm™ heating/stirring module (uniform dry heat to the sample) coupled with the Thermo Scientific™ Reacti-Vap™ evaporator (simultaneous or separate delivery of pressurized gas) to provide a complete solution for derivatization or other small scale reactions.

Reacti-Therm Tips & Tricks

The Reacti-Therm temperature is set by turning to heat knob until the desired temperature appears in the LED display

Tip 1

The displayed value will continue to flash until the Reacti-Block temperature approaches the set temperature. Allow 30-60 minutes to reach the set temperature, then a further 15 minutes for the temperature to stabilise



Reacti-Therm units are available in single and triple block configurations to hold 1 or 3 Reacti-Blocks respectively.

Tip 2

If using the triple block version, the unit must be filled with three Reacti-Blocks in order to obtain acceptable heating performance



Reacti-Therm units are available with heating only and heating / stirring functions

Tip 3

Select the heating / stirring functions if working with samples that are viscous or sticky in nature. Here, a magnetic stirrer bar needs to be inserted into the vial



Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand

A more accurate and responsive control of the heating function of the Reacti-Therm can be achieved

Tip 4

Use a Remote Temperature Probe in conjunction with a standard thermometer

The maximum pressure specification for a Reacti-Vap is 2PSi

Tip 5

Each Reacti-Vap has a test run for a maximum pressure, but not a maximum flow rate. Reacti-Vaps are pressure rated to a maximum of 2PSI. There is a pressure relief valve connected to the side of the Reacti-Vap that will “pop off” when the internal pressure reaches 2PSI, not allowing it to go any higher.

Reacti-Vap parts can easily be replaced

Tip 6

A range of replacement parts, such as needles, luer-lok® plugs and mounting brackets are available and can easily be changed by a user



Reacti-Therm Dry Block Sample Incubation Systems perform derivatization and other small-scale reactions reliably and with ease. It features Reacti-Therm Modules that utilize dry block heaters to provide uniform, stable heating with a temperature range of 10–200 °C. The Reacti-Therm Dry Block Sample Incubation System is available in four models—single block and triple block sizes with either heat-only or heat-and-stir capability. This versatile system is compatible with a range of interchangeable accessories and is ideal for performing heating, stirring, and evaporation functions on small samples for a range of applications.

See more at: www.thermoscientific.com/chromatography

For more information, visit our website at www.thermoscientific.com/chromatography

© 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

USA and Canada +1 800 332 3331
Australia 1300 735 292 (free call domestic)
China 800 810 5118 (free call domestic)
400 650 5118
France +33 (0)1 60 92 48 34
Germany +49 (0) 2423 9431 20 or 21
India +91 22 6742 9494
+91 27 1766 2352

Japan 0120 753 670 (free call domestic)
0120 753 671 (fax)
Korea +82 2 3420 8600
United Kingdom +44 (0) 1928 534 110
New Zealand 0800 933 966 (free call domestic)
Singapore +65 6289 1190
All Other Enquiries +44 (0) 1928 534 050

Technical Support

For advice and support, please visit our website:
www.thermoscientific.com/chromexpert

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand