

Determinations of Targeted Metabolites Using Capillary Ion Chromatography High Resolution Accurate Mass Spectrometry

Terri Christison,¹ Junhua Wang,² Linda Lopez,¹ Yingying Huang²

¹Thermo Fisher Scientific, Sunnyvale, CA, USA; ²Thermo Fisher Scientific, San Jose, CA, USA



Overview

Purpose: To develop and evaluate a quantitative method for targeted metabolic profiling analysis of mono-, di-, and tri-phosphate nucleotides using capillary ion chromatography tandem mass spectrometry (Cap IC-MS-MS).

Methods: A 5 μ L sample was separated on a Thermo Scientific™ Dionex™ IonSwift™ MAX-100 monolith anion-exchange capillary column using an hydroxide gradient elution. Target analytes were detected using a Thermo Scientific™ Quantum TSQ Access MAX™ Triple Quadrupole Mass Spectrometer in selected reaction monitoring (SRM) mode. Potassium hydroxide mobile phase was converted to water, post-column by a Thermo Scientific™ Dionex™ ACES™ 300 Anion Capillary Electrolytic Suppressor. Acetonitrile was added to the IC stream to assist in desolvation and to improve sensitivity.

Results: Target nucleotides were chromatographically resolved within 40 min; two SRM transitions were used for quantitation and confirmation for each analyte. The chromatographic separation was essential to eliminate the SRM interferences from structurally related analytes, e.g. ADP and ATP. The Limit of Quantitation was 1 nM for each analyte (5 fmol on column). Calibration range was 1–1000 nM with a coefficient of determination (r^2) > 0.99. Good quantitation accuracy and precision at 10 nM, 100 nM and 500 nM, were also achieved, 86.4% to 107% and %RSD less than 6% (< 17.4% for 10 nM), respectively.

Introduction

Nucleotides are essential compounds active in many cell functions such as energy storage and release, signal transduction and synthesis of DNA and RNA. In recent years, there have been extensive studies of using nucleoside analogs as drugs in anti-cancer, anti-viral and immunosuppressive therapy,^{1,2} and monitoring of their activated nucleotides metabolites is of paramount importance to understand the pharmacology. Reported methods for nucleotides quantification include liquid chromatography (LC) with ion pairing reagent,³ hydrophilic interaction liquid chromatography⁴ (HILIC), ion chromatography^{5,6} (IC) and capillary electrophoresis^{7–9} (CE) with different detections such as conductivity, UV and mass spectrometry (MS). These methods usually do not have sufficient retention or complete chromatographic resolution to establish nucleotides profile, or required sensitivity to quantify nucleotides at trace levels.

This study describes a capillary IC tandem MS method for profiling analysis of 19 native and two modified nucleotides. Chromatographic separation was achieved on a monolith column with baseline resolutions for most of the target analytes. MS was operated in selected reaction monitoring (SRM) mode to achieve sensitive and selective quantitation. Ion exchange column chemistry with the monolith format provides the required chromatographic selectivity and the speed for the analyses. The reduced flow rate used in the method further improves the MS detection sensitivity, thus ensuring the accurate quantitation at low nM level. Three isotope labeled internal standards were used for better quantitation accuracy. Potassium hydroxide mobile phase was converted to water by a suppressor ensuring compatibility with MS detection.

Methods

Sample Preparation

HEK 293T Cells were cultured in ATCC-recommended medium at 37 °C and in 5% CO₂ atmosphere. All media were supplemented with 10% fetal bovine serum, 100 I.U./mL penicillin, and 100 μ g/mL streptomycin. After growing to 70% confluence (at a density of 106 cells/mL), cells were washed with PBS twice then cultured in PBS containing 0 μ M and 250 μ M methylglyoxal for 3 hours.

After methylglyoxal treatment, cell pellets were harvested from the drug- or mock-treated cells (~ 2 \times 10⁷ cells), washed with PBS buffer and re-suspended in 10 mM sodium citrate (pH 4.5) at a final volume of 400 μ L.

The mixture was sonicated for 1 min, to which solution was added 4 μ L of 50 mM DTT and 50 μ L 100 mM D-P. The resulting solution was incubated at 37 °C for 30 min. After centrifugation at 13,000 rpm for 5 min, the supernatant was transferred out and filtrated using 3000-MW cut-off Microcon® ultra centrifugation units. The filtrates were then subjected to IC-MS-MS analysis.


Ion Chromatography

- Thermo Scientific™ Dionex™ ICS-5000+ Capillary HPIC™ system (capillary DP Dual Pump, EG Eluent Generator, DC Detector Column modules)
- Thermo Scientific™ Dionex™ IC Cube™ with 6-port injection valve
- Thermo Scientific Dionex AS-AP Autosampler

Columns: Dionex IonSwift MAX 100 with guard, 0.25 × 250 mm
Eluent Source: Thermo Scientific Dionex EGC-KOH capillary cartridge
KOH Gradient: 10 mM (-7–0 min), 10–75 mM (0.1–25 min),
75 mM (25–29.9 min), 75–100 mM (29.9–30 min), 100 mM (30–39.9 min),
100–10 mM (39.9–40 min)
Flow Rate: 0.015 mL/min
Inj. Volume: 5 μL
Column Temp.: 40 °C

All tubing were precision cut to avoid extra dead volumes. A MicroTee® (Upchurch Scientific, P-890) was used to blend the IC stream with acetonitrile desolvation solvent before entering ESI ionization interface. Chromatography is not shown. The preferred Cap IC-MS-MS flow schematic is shown in Figure 1.

FIGURE 1. Preferred Cap IC-MS-MS schematics.



Mass Spectrometry

A Thermo Scientific™ Quantum TSQ Access MAX™ Triple Quadrupole Mass Spectrometer was used in this study and coupled to the Capillary RFIC IC system with a Ion Max source and heated ESI probe (HESI II) with low-flow metal capillary kit. The source parameters were: spray voltage (3500 V), vaporizer temperature (150 °C), sheath gas pressure (25 arbitrary unit), aux gas pressure (15 arbitrary unit), capillary temperature (200 °C). Two SRM transitions were used for the quantitation (Q-SRM) and confirmation (C-SRM) of each target analyte with collision energy (CID) optimized for each SRM transition. Detailed SRM scan events are listed in Table 1.

Data Analysis

Thermo Scientific™ Xcalibur™ version 2.2 SP1 with Foundation version 2.0 SP1 and TSQ version 2.3 SP3.

Thermo Scientific™ Dionex™ DCMSLink™ version 2.11.

Results

Method Performance

As seen in Figure 2 of the SRM chromatograms, all target analytes were separated with necessary chromatographic resolution, and closely eluted peaks were successfully differentiated by SRM.

As summarized in Table 2, the calibration for each analyte from 1 nM to 1000 nM were very good, $r^2 > 0.99$. The accuracy and precision were evaluated at 10, 100, and 500 nM levels. Accuracy was observed in the range from 86.4% (GMP, 10 nM) to 107% (GDP, 100 nM). Precision was addressed as %RSD and observed in the range from 1.65% (IDP, 100 nM) to 17.4% (UDP, 10 nM).

TABLE 1. SRM Scan Events and Quantified Nucleotides in Cell Extracts

R.T. (min)	Analyte	Precursor	Q-SRM	C-SRM	Treated*	Untreated*
3.2	CMP	322.0	79 (40)	97 (24)	> 40	> 40
5.5	AMP	346.0	79 (37)	134 (36)	> 100	> 100
10.7	CDP	402.0	159 (25)	384 (20)	5	8.66
13.0	UMP	323.0	79 (37)	97 (25)	> 100	> 100
15.3	ADP	426.0	328 (19)	159 (27)	13.8	21.2
17.2	IMP	347.0	79 (36)	135 (32)	37.2	23.8
17.4	GMP	362.0	211 (21)	79 (37)	> 100	> 100
18.0	dCTP	466.0	159 (32)	368 (23)	ND	ND
19.4	CTP-d ₉	491.0	159 (28)	393 (19)		
19.5	CTP	482.0	159 (27)	384 (21)	0.351	1.17
22.8	UDP	403.0	159 (27)	111 (22)	7.64	11.6
23.4	dATP	490.0	159 (26)	392 (25)	ND	ND
24.4	ATP-d ₁₀	516.0	159 (31)	418 (25)		
24.4	ATP	506.0	408 (23)	159 (30)	0.524	2.53
27.6	IDP	427.0	329 (19)	135 (25)	ND	ND
27.8	GDP	442.0	344 (20)	150 (27)	9.05	10.1
29.2	dTTP	481.0	159 (27)	383 (20)	0.085	0.100
30.7	UTP-d ₉	492.0	159 (31)	394 (17)		
30.7	UTP	483.0	385 (22)	159 (36)	0.227	1.52
33.2	dGTP	506.0	159 (36)	408 (22)	ND	ND
33.7	ITP	507.0	409 (21)	159 (36)	ND	ND
33.8	GTP-d ₁₀	532.0	159 (27)	434 (22)		
33.8	GTP	522.0	424 (22)	159 (27)	0.374	1.45
36.1	N ₂ -CEdGTP	578.0	480 (22)	159 (40)	ND	ND
37.0	N ₂ -CEGTP	594.0	496 (22)	159 (40)	ND	ND

FIGURE 2. Profiling 19 native and two modified nucleotides using Cap IC-MS-MS.




TABLE 2. Method Performance

Analyte	r ²	10 nM			100 nM			500 nM		
		Ave	%RSD	%Acc*	Ave	%RSD	%Accuracy	Ave	%RSD	%Acc*
CMP	0.998	9.91	8.6	99.1	98.5	4.2	98.5	476	3.6	95.2
AMP	0.998	9.72	10	97.2	101	3.8	101	507	1.8	101
CDP	0.994	9.45	6.3	94.5	94.0	3.2	94.0	478	2.1	95.6
UMP	0.997	9.70	9.9	97.0	97.9	2.9	97.9	526	3.7	105
ADP	0.995	10.2	8.4	102	96.7	4.0	96.7	504	2.3	101
GMP	0.995	8.64	15	86.4	106	4.4	106	519	2.4	104
dCTP	0.999	10.2	4.6	102	98.2	2.8	98.2	497	2.1	99.4
IMP	0.995	8.77	15	87.7	92.7	3.5	92.7	472	4.6	94.4
CTP	0.999	10.4	11	104	97.0	3.1	97.0	487	1.5	97.4
UDP	0.997	10.1	17	101	99.2	4.6	99.2	533	3.8	107
dATP	0.999	10.5	8.9	105	96.8	1.8	96.8	510	1.7	102
ATP	0.998	10.1	16	101	96.4	2.6	96.4	500	1.8	100
GDP	0.994	10.2	12	102	104	1.8	104	534	2.7	107
dTTP	0.998	9.95	11	99.5	99.6	3.0	99.6	544	3.3	109
IDP	0.997	9.48	12	94.8	99.0	1.7	99.0	522	3.5	104
UTP	0.998	9.62	16	96.2	94.8	2.3	94.8	519	2.6	104
dGTP	0.997	10.5	9.9	105	98.8	3.0	98.8	516	2.3	103
GTP	0.996	10.5	12	105	97.5	4.14	97.5	502	4.0	100
ITP	0.997	9.96	9.4	99.6	98.1	5.8	98.1	498	3.8	99.6

*Acc = % Accuracy

The treated and untreated cell extracts were diluted 100-fold with deionized water and 5 μ L of diluted samples were injected for Cap IC-MS-MS analysis. Most native nucleotides were detected (Table 1). The mono-phosphate nucleotides were the most abundant ones. However, modified triphosphate nucleotides were not detected in either sample probably due to degradation to di- or mono-phosphates.

Conclusion

A Cap IC-MS-MS method was developed for the targeted profiling quantitation of nucleotides in biological samples. This method was evaluated and demonstrated:

- Excellent chromatographic separation for target analytes
- Selective SRM MS-MS detection capable of differentiating closely eluted analytes
- Ultra-sensitive detection of target analytes down to 1 nM with only few μ L sample consumption
- Successful application for biological sample analysis.

References

1. Cohen, S.; Jordheim, L. P.; Megherbi, M.; et al. *J. Chromatogr. B* **2010**, 878 (22), 1912-1928.
2. Jansen, R. S.; Rosing, H.; Schellens, J. H. M.; Beijnen, J. H. *Mass Spectrom. Rev.* **2011**, 30 (2), 321-343.
3. Seifar, R. M.; Ras, C.; van Dam, J. C.; et al. *Anal. Biochem.* **2009**, 388 (2), 213-219.
4. Johnsen, E.; Wilson, S. R.; Odsbu, I.; Krapp, A.; et al. *J. Chromatogr. A* **2011**, 1218 (35), 5981-5986
5. C. Caldwell, I. *J. Chromatogr. A* **1969**, 44 (0), 331-341.
6. Inoue, K.; Obara, R.; Akiba, T.; Hino, T.; et al. *J. Agric. Food. Chem.* **2008**, 56 (16), 6863-6867
7. Nguyen, A. L.; Luong, J. H. T.; Masson, C. *Anal. Chem.* **1990**, 62 (22), 2490-2493.
8. Soga, T.; Ohashi, Y.; Ueno, Y.; Naraoka, H.; et al. *Journal of Proteome Research* **2003**, 2 (5), 488-494

Acknowledgements

The HEK 293T cells were provided from Yinsheng Wang as part of his doctoral research at the University of California, Riverside, CA, USA

www.thermoscientific.com

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. Microcon is a registered trademark of EMD Millipore, and MicroTee is a registered trademark of Upchurch. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Fisher Scientific, Sunnyvale, CA
USA is ISO 9001:2008 Certified.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Brazil +55 11 3731 5140
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 6 6885 1213
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
Taiwan +886 2 8751 6655
UK/Ireland +44 1442 233555
USA +1 800 532 4752

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific