



Thermo SCIENTIFIC

Potentiometric Titration Application Notes

Applications Log # 543

Overview

The concentrations of sodium and % sodium chloride were determined by a Thermo technique called KAP Analysis. Aliquots of sodium standard are added automatically to a sample containing an Orion sodium electrode. The Orion 960 Autochemistry System calibrates the electrode and calculates the concentration.

Market	Food and Beverage	Species Measured	Sodium
Sample	Meat Flavorings	Sample Size	1.0g
		Typical Concentration	0.05% wt/wt
Technique #	2 Multiple Known Addition	Electrode	Sodium 8611BN
Solutions	Sodium 1M KAP 650700; Sodium ISA 841111. Electrode fill solution 900010		
Sample Prep	Weigh out 20g sample & transfer to 100ml volumetric flask. fill w/ DI H ₂ O. Pipette 5ml aliquot of sample place in a 100ml beaker, add 50ml DI water and 5ml sodium ISA. Titrate.		
Statistics			
# of Trials	4	Mean	0.0524%wt/wt
		%CV	1.31
		Analysis Time	2.0minute(s)
Comments	Rinse the electrodes, stirrer, and dispenser probe between measurements with deionized water.		

Method Parameters

Sample Volume/Weight	1.00 g	Timed or Stability Readings	18.0 mV/min stability
Constant Increment	0.196 mL	Number of Endpoints	1
Max Titrant Volume	10.00 mL	Desired Units	% w/w
Molecular weight	22.99 g for beef c anion	Predose	
Prestir	120.0 second(s)	Additional Parameters	
Reaction Ratio	1.00		