

Robust and Sensitive Relative Quantification of Isobarically Labeled Glycans by LTQ Orbitrap Velos Hybrid Mass Spectrometer

Rosa Viner¹, David M. Horn¹, Shujuan Tao², Ron Orlando²

¹Thermo Fisher Scientific, San Jose, CA, ²Complex Carbohydrate Research Center, University of Georgia, Athens, GA, USA

Overview

Purpose: Evaluate the performance of a new hybrid ion trap – Orbitrap mass spectrometer featuring a dual-pressure ion trap and an improved HCD collision cell for MS and MSⁿ-based QUIBL quantification of *N*-linked glycans.

Methods: The *N*-linked glycans were permethylated with ¹³CH₃I or ¹²CH₂DI, mixed in different ratios, analyzed by MS/MSⁿ on a LTQ Orbitrap Velos, and quantified.

Results: Isobarically labeled glycans mixed in different ratios were successfully quantified using full MS, HCD or CID/HCD MSⁿ spectra over 2 orders of magnitude. The maximum resolution of the Orbitrap mass spectrometer, 100,000 at *m/z* 400, was required to resolve and quantify the ¹³CH₃ and ¹²CH₂D species.

Introduction

A current challenge in the field of glycomics is determining how to quantify changes in glycan expression between different cells, tissues, or biological fluids. Atwood *et al*¹ described a novel strategy for relative quantification of glycans by isobaric labeling (QUIBL) that requires high resolution and mass accuracy to distinguish isobaric pairs. In this report, the performance of a new hybrid ion trap – Orbitrap mass spectrometer featuring a dual-pressure ion trap and an improved HCD collision cell² is evaluated for MS and MSⁿ-based QUIBL quantification of fetuin *N*-linked glycans and lacto-*N*-fucopentaoses isomers.

Methods

Sample Preparation

Bovine fetuin *N*-linked glycans were released by overnight incubation with PNGase F (Sigma), separated from peptides by reverse phase chromatography, permethylated with ¹³CH₃I or ¹²CH₂DI (Cambridge Isotopes), and then dried and redissolved in 50% MeOH and 1mM NaOH for MS analysis¹. Isobarically labeled Lacto-*N*-fucopentaose I (LNFP I), Lacto-*N*-fucopentaose III (LNFP III), Lacto-*N*-difucohexaose I (LNDFH I), and Lacto-*N*-difucohexaose II (LNDFH II) from V-labs were prepared as described for the fetuin glycans. Mixtures of glycans were prepared with concentration ratios ranging from 20:1 to 1:10.

MS Analysis

Glycans at pmol concentrations were analyzed on a Thermo Scientific LTQ Orbitrap Velos hybrid mass spectrometer via direct infusion using static nanospray. Data were acquired in full MS, HCD, and FT CID MSⁿ scan modes at 60,000-100,000 resolving power (RP) at *m/z* 400.


Data Processing

Thermo Scientific Xtract software was used for spectral deconvolution and monoisotopic mass determination. Quantification was performed by summing the ¹³CH₃ or ¹²CH₂D-labeled total ion intensities for each glycan and applying correction factors for permethylation completeness.

Results

Figure 1 illustrates the QUIBL method. Bovine fetuin *N*-linked glycans are permethylated with either ¹²CH₂DI or ¹³CH₃I (Figure 1a, 1b) to generate isobaric pairs. The permethylated samples were mixed in different ratios and analyzed using an LTQ Orbitrap VelosTM mass spectrometer, as high resolution mass spectrometer is required. The mass difference between the isobaric labels is small (0.00292 Da) but glycans contain multiple permethylation sites (shown in parentheses in Figure 2) that together produce a mass difference that is completely resolved using 100,000 RP, but not 60,000 at *m/z* 400 (Figure 1d versus 1c). Relative quantification by QUIBL can be performed either in MS mode (Figure 1) or by MS/MS performed at high resolution for analysis of isomeric glycan mixtures (Figure 5b). For accurate relative quantification by QUIBL, correction factors for each label need to be calculated. As shown in Figure 1a versus 1b, the ¹²CH₂D labeled sample shows lower level isotopic of isotopic enrichment than the ¹³CH₃ labeled sample (97.6% versus 98.9%) plus more incomplete permethylation due to a kinetic isotope effect³.

FIGURE 1. The QUIBL method is demonstrated for the fetuin tri-antennary glycan permethylated with ¹²CH₂DI (a) or ¹³CH₃I (b), and its 1:1 mixture. FTMS spectra were acquired at 60K (a, b, c) or 100K (d).



Nevertheless, it is still possible to perform QUIBL quantification using the sum of all fragment isotope ion abundances as described in *Methods*. For the 1:1 mixture for fetuin triantennary glycan at 100,000 RP (Figure 1d), a ratio of 1.16 was observed.

The dynamic range and linearity of QUIBL on the LTQ Orbitrap Velos mass spectrometer were evaluated using 4 main species of bovine fetuin *N*-glycans (Figure 2), mixed in 5 different ratios (20:1-1:10) using 100K resolving power in the MS and MS/MS modes using HCD as shown in Figures 3 and 4. These results demonstrate that quantification using the QUIBL method is linear over at least 2 orders of magnitude, even with incomplete isotopic incorporation. Sialic acid ions generated by HCD provided the most accurate results for quantification for the higher ratios due to highest resolution achieved (>100 K, Figure 4). The overall maximum error was below 15%, which is comparable to other isotopic based quantification techniques.

FIGURE 2. Orbitrap MS spectrum (a) and HCD spectra (b) of bovine fetuin *N*-linked permethylated glycans mixture and its compositions acquired at 100K RP at *m/z* 400.




FIGURE 3. FTMS QUIBL quantification of the fetuin bi-antennary glycan. Orbitrap spectra were acquired at 100K RP. Inset shows correlation between the theoretical and experimental ratios for mixtures at 20:1, 10:1, 1:1, 0.5:1, and 1:10 (¹³CH₃:CH₂D).






FIGURE 4. HCD QUIBL quantification of the fetuin tetra-sialylated glycan using sialic acid fragment ion at *m/z* 382. Orbitrap spectra were acquired at 100K RP. Inset shows correlation between the theoretical and experimental ratios for mixtures at 20:1, 10:1, 1:1, and 0.5:1 (¹³CH₃:CH₂D).



One of the unique aspects of QUIBL is that simultaneous quantification of isomeric glycans can be performed using MSⁿ mode. Two pairs of glycan isomers, LNFP I and LNFP III, LNDFH I and LNDFH II, were used to demonstrate such quantification. Each glycan was permethylated with ¹³CH₃I or ¹²CH₂DI and mixed in different ratios. Results for the LNFP I and LNFP III pair mixed in 1:5:5:1 ratio are shown in Figure 5. Isotope enrichment for these experiments was calculated as 98.9% for ¹³CH₃ and 98.1% for ¹²CH₂D. LNFP I and LNFP III have exactly the same molecular mass ([M+2Na]²⁺ = 569.8) and, as expected, the ratio calculated using the full MS spectra was close to 1 (Figure 5a).

An isomeric mixture such as LNFP I and LNFP III requires MS/MS performed at high resolution for correct quantification. Introduction of an integrated C-trap/HCD collision cell combination² for the LTQ Velos Orbitrap mass spectrometer significantly increases the sensitivity of the instrument in HCD mode, facilitating QUIBL analysis. The ion at *m/z* 569.8 ([M+2Na]²⁺) was selected and fragmented by HCD. The ion at *m/z* 439.2 was used for relative quantification of LNFP I. The signature ion at *m/z* 907.4 was used for quantification of LNFP III as shown in Figure 5b (Insets). Both calculated glycan ratios were within 10% of their theoretical values. LNDFH I and LNDFH II demonstrated the same relative performance (data not shown). Overall similar dynamic range and linearity of the QUIBL method was obtained for isomer mixtures quantification as for fetuin *N*-glycans.

FIGURE 5. Orbitrap MS (a) and HCD (b) spectra of isomeric mixture of isobarically labeled LNFP I and LNFP III in 1:5:5:1 ratio. Orbitrap spectra were acquired at 100K RP. Insets show zoom in of ions used for quantification.



Conclusion

- The QUIBL method was successfully used for relative quantification of glycans in MS or MS/MS mode on a LTQ Orbitrap Velos mass spectrometer over 2 orders of magnitude within <15% error.
- Optimization of the reaction conditions for permethylation and correction factors for rate of incorporation ¹³CH₃I versus ¹²CH₂DI is required.
- Successful quantification of isomer mixtures by HCD at 100,000 FWHM was demonstrated.
- QUIBL was validated as a powerful quantitative glycomics approach.

References

- Atwood, J.; Cheng, L.; Alvarez-Manilla, G.; *et al.*, (2007) Quantitation by isobaric labeling: applications to glycomics. *Journal of Proteome Research* 7, 367-374.
- Olsen, J.; Schwartz, J.C.; Griep-Raming, J.; *et al.*, (2009) A dual pressure linear ion trap – Orbitrap instrument with very high sequencing speed. *Mol. Cell. Proteomics* 8, 2759-2769.
- Botelho, J.; Atwood, J.; Cheng, L.; *et al.*, (2008) Quantification by isobaric labeling (QUIBL) for the comparative glycomics study of O-linked glycans. *Int. Journal of Mass Spectrometry* 278, 137-142.

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.