



Salmonella Rapid Culture Method

Confidence in food safety:
2 steps, 2 days

[Contact Me](#)

[Rapid vs Traditional](#)

[Product Information](#)

[AOAC License](#)

[AOAC Journal Article \(pdf\)](#)

[AOAC Annual Conference](#)

[Presentation](#)

[Press Release](#)

[Remel Home](#)

The Salmonella Rapid Culture Method

- Simple two-step procedure — one broth to one plate
 - **Step 1** Enrichment: ONE Broth-Salmonella
 - **Step 2** Isolation: Brilliance™ Salmonella
- Reduced time to results — less than two days
- No broth agitation required
- Chromogens in Brilliance™ Salmonella aid in easy identification and differentiation
- *Salmonella* spp. are visible as purple colonies

Rely on rapid results from Remel for successful Salmonella detection!

[Clinical Products](#) | [Industrial Products](#) | [Support](#) | [Announcements](#) | [About Remel](#) | [Careers at Remel](#) | [Contact Us](#) | [International](#) | [Site Map](#) | [Terms of Use and Privacy Notice](#) | [Terms and Conditions of Sale](#) | Thermo Fisher Scientific Inc.

12076 Santa Fe Trail Drive, Lenexa, KS 66215 USA | 800.255.6730
©2009 Thermo Fisher Scientific. All Rights Reserved